These guidelines are for IMS internal use only and are not for public dissemination. IMS competency validation guides are strictly advisory and confer no certifications or warranties. These evaluations are intended to train IMS sales representatives to focus on personnel, processes, environments and equipment that may need attention. These competency guides are based on original equipment manufacturer information and related professional society standards for practice and are not a substitute for facility personnel in developing and/or meeting the requirements of accrediting or regulatory elements of performance.

Mechanical Leak Test Competency Verification Checklist



Employee Name / Title Territory/Region Validator's Name			A Surgical Device Management & Clinical Consulting Company Integrated Medical Systems International, Inc. www.imsready.com: 800.783.9251				
					•	Needs Assistance	Correction Plan
				Removes detachable parts of the scope.			
Checks light guide prong and housing for tightness.		o l					
Turns on maintenance unit.		o l					
Fills sink with water to completely cover.		0					
Attaches water resistant cap to scope.		o l					
Attaches leak tester to scope.							
Insufflates scope prior to submersion assuring that bending rubber has expanded.	□						
Flushes water though all scope channels with syringe to remove trapped air.							
Observes scope for at least 90 seconds.							
Angulates scope in all directions while checking bending rubber for leaks.	□						
Angulates scope in all directions while checking biopsy channel for leaks.	٥						
Angulates scope in all directions while checking control knobs for leaks.	□						
Angulates elevator raiser evaluating lifter and elevator		0					
Inspects for leaks – fluid resistant cap.							
Inspects for leaks – insertion tube.		0					
Inspects for leaks – light guide connector and tube.		0					
Inspects for leaks - control head.							
Massages all video switches with a circular motion.		0					
Keeps scope completely submerged for entire leak test.							
Removes scope from sink.		0					
Turn off maintenance unit.							
Unplug leak tester.		0					
Disconnect leak tester from scope.							
Wipes excess liquid from scope.		0					
Can access procedure for further processing a scope with an identified leak.	┚						
Action Plan:							
Performer's Signature			Todav's Date				

Validator's Signature____

Expiration Date